

Anti-IL5RA antibody



Catalog Number: 141417

Product name

Anti-IL5RA antibody

Immunogen

[Mouse IL5RA \(CD125\) \(His Tag\) recombinant protein](#)

Specificity

Mouse IL5RA (CD125)

Antibody description

Rabbit Polyclonal to Mouse IL5RA

Preparation

Produced in rabbits immunized with purified, recombinant Mouse IL5RA (CD125). IL5RA (CD125) specific IgG was purified by Mouse IL5RA (CD125) affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 5-20 µg/ml

Validations



Lanes	A
Items	
Sample	Mouse heart tissue lysate
Sample Volume (µg/lane)	30
Gel	13% SDS-PAGE reducing gel
Recommended Concentration	5-20 µg/ml
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.
Developed using Odyssey imaging system.	
Explanation	Predicted band size : 48 kDa Observed band size : 44 kDa

Mouse IL5RA (CD125) Western blot (WB) 14904