

# PGC1 alpha Polyclonal Antibody



Catalog Number: 163417

## Product name

Rabbit IgG

PGC1 alpha Polyclonal Antibody

## Applications

## Specificity

WB

Human, Mouse, Rat

## Dilutions

## Antibody description


WB 1:500 - 1:2000

Polyclonal antibody to PGC1 alpha

## Validations

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of mouse PGC1 alpha (NP\_032930.1).



## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Western blot - PGC1 alpha Polyclonal Antibody

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of U-87MG cells, using PGC1 alpha antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

## Clonality

Polyclonal

## Ig Type