

General Information

Protein Construction

A DNA sequence encoding the Influenza A virus (A/duck/NZL/76/1984(H9N1)) hemagglutinin (Q20PC1) (Met1-Arg388), termed as HA1, was expressed with a polyhistidine tag at the C-terminus.

Organism

H9N1

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Asp 19

Molecular Mass

The recombinant HA1 subunit of the Influenza A virus (A/duck/NZL/76/1984(H9N1)) hemagglutinin consists of 331 amino acids and predicts a molecular mass of 37 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.