

# GAD2 antibody

Catalog Number: 110817

## Product name

GAD2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GAD2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Positive IP detected in rat brain tissue. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of GAD2 recombinant protein (Accession Number: NM\_001134366). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110817(GAD2 Antibody) at dilution of 1:600



IP Result of anti-GAD2 (IP:Catalog No:110817, 3ug; Detection:Catalog No:110817 1:300) with rat brain tissue lysate 4000ug.