

DUSP15 antibody



Catalog Number: 110117

Product name

DUSP15 antibody

Immunogen

[Human DUSP15 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

DUSP15 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Observed molecular weight by Western-blot: 32kd

Preparation

This antibody was obtained by immunization of DUSP15 recombinant protein (Accession Number: NM_080611). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110117(DUSP15 antibody) at dilution of 1:300