

KPTN antibody



Catalog Number: 112117

Product name

KPTN antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity


Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

KPTN Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human cerebellum tissue, mouse skeletal muscle tissue, NIH/3T3 cells. Positive IHC detected in human kidney tissue, human heart tissue, human lung tissue, human ovary tissue, human placenta tissue, human skin tissue, human spleen tissue. Observed molecular weight by Western-blot: 48 kDa

Validations




Preparation

This antibody was obtained by immunization of KPTN recombinant protein (Accession Number: NM_007059). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112117(KPTN antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Immunohistochemical of paraffin-embedded human kidney using Catalog No:112117(KPTN antibody) at dilution of 1:100 (under 10x lens)



Dilutions

Immunohistochemical of paraffin-embedded

KPTN antibody



Catalog Number: 112117

human kidney using Catalog No:112117(KPTN antibody) at dilution of 1:100 (under 40x lens)