

RLN2 Polyclonal Antibody



Catalog Number: 165117

Product name

RLN2 Polyclonal Antibody

Specificity

Human, Rat

Antibody description

Polyclonal antibody to RLN2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-185 of human RLN2 (NP_604390.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations



Western blot - RLN2 Polyclonal Antibody

Western blot analysis of extracts of 22Rv1 cells, using RLN2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time:



Immunohistochemistry - RLN2 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using RLN2 antibody at dilution of 1:200 (40x lens).