

# Anti-FOXF2 antibody



Catalog Number: 175117

## Product name

Anti-FOXF2 antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to FOXF2

## Dilutions

WB: 1:500

## Preparation

This antigen of this antibody was peptide

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 50% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Validations




Fig1: Western blot analysis of FOXF2 on human placenta tissue lysate using anti-FOXF2 antibody at 1/200 dilution.