

Anti-FOXF2 antibody



Catalog Number: 175117

Product name

Anti-FOXF2 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to FOXF2

Preparation

This antigen of this antibody was peptide

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 50% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500

Validations




Fig1: Western blot analysis of FOXF2 on human placenta tissue lysate using anti-FOXF2 antibody at 1/200 dilution.