

Anti-EpCAM Mouse antibody



Catalog Number: 169117

Product name

Anti-EpCAM Mouse antibody

Specificity

Human

Antibody description

EpCAM Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human EPCAM (AA: Extra(24-265)) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

FC

Validations




Western blot analysis using EPCAM mAb against



human EPCAM (AA)

Western blot analysis using EPCAM mouse mAb against A431 (1), MCF-7 (2) cell lysate.




Western blot analysis using EPCAM mAb against HEK293 (1) and EPCAM (AA)

Anti-EpCAM Mouse antibody



Catalog Number: 169117



Flow cytometric analysis of A431 cells using EPCAM mouse mAb (green) and negative control (red).