

# Anti-EpCAM Mouse antibody



Catalog Number: 169117

## Product name

Anti-EpCAM Mouse antibody

## Specificity

Human

## Antibody description

EpCAM Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human EPCAM (AA: Extra(24-265)) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

FC

## Validations




Western blot analysis using EPCAM mAb against



human EPCAM (AA)

Western blot analysis using EPCAM mouse mAb against A431 (1), MCF-7 (2) cell lysate.




Western blot analysis using EPCAM mAb against HEK293 (1) and EPCAM (AA)

# Anti-EpCAM Mouse antibody



Catalog Number: 169117

---



Flow cytometric analysis of A431 cells using EPCAM mouse mAb (green) and negative control (red).