

CDCA7 antibody



Catalog Number: 109117

Product name

CDCA7 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CDCA7 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 40-43 kDa

Preparation

This antibody was obtained by immunization of CDCA7 recombinant protein (Accession Number: NM_145810). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:10-1:100

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109117(CDCA7 Antibody) at dilution of 1:200



Immunofluorescent analysis of HepG2 cells using Catalog No:109117(CDCA7 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG