

# Mouse TNFSF12 (Fc Tag) recombinant protein



Catalog Number: 500517

## General Information

### Gene Name Synonym

TNF-related weak inducer of apoptosis; Tumor necrosis factor ligand superfamily member 12, membrane form; Tumor necrosis factor ligand superfamily member 12, secreted form

### Protein Construction

A DNA sequence encoding the mouse Tnfsf12 (NP\_035744.1) (Arg105-His249) was expressed with the Fc region of human IgG1 at the N-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

> (83.8+3.6) % as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Glu

## Molecular Mass

The recombinant mouse Tnfsf12 consists 405 amino acids and predicts a molecular mass of 44.4 kDa.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Tnfsf12 protein SDS-PAGE