

# Anti-HA antibody

Catalog Number: 142517

## Product name

Anti-HA antibody

## Immunogen

[H1N1 HA \(His Tag\) recombinant protein](#)

## Specificity

Influenza A H1N1 (A/WSN/33) Hemagglutinin / HA1

## Antibody description

Rabbit Polyclonal to H1N1 HA

## Preparation

Produced in rabbits immunized with purified, recombinant Influenza A H1N1 (A/WSN/33) Hemagglutinin / HA1 . Influenza A H1N1 (A/WSN/33) Hemagglutinin / HA1 specific IgG was purified by Influenza A H1N1 (A/WSN/33) Hemagglutinin / HA1 affinity chromatography.

## Formulation

0.2 µm filtered solution in PBS

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of

properly. Avoid repeated freeze-thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, ELISA, ICC/IF, IF, IP


## Dilutions

WB: 0.5-5 µg/ml

ELISA: 0.1-0.2 µg/ml

This antibody can be used at 0.1-0.2 µg/ml with the appropriate secondary reagents to detect Influenza A H1N1 (A/WSN/33) Hemagglutinin / HA1.

## Validations



Items	Lanes	A	B
Sample (Recombinant protein)	H1N1-WS-33-HA1-His (Cat#11692-V08H1)		
Sample Volume (ng/lane)		50	10
Gel	13% SDS-PAGE reducing gel		
Recommended Concentration	0.5-5 µg/ml		
Secondary Antibody	Dylight 800-labeled Antibody To Rabbit IgG (H+L), at 1:5000 dilution.		
Developed using Odyssey imaging system.			

Influenza A H1N1 (A/WSN/33) Hemagglutinin / HA1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot