

MT2A Polyclonal Antibody



Catalog Number: 163517

Product name

MT2A Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to MT2A

Preparation

Antigen: Recombinant protein of human MT2A.

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

IF

Dilutions

IF 1:50 - 1:200

Validations

Immunofluorescence - MT2A Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using MT2A antibody at dilution of 1:100. Blue: DAPI for nuclear staining.