

# Rad51D antibody



Catalog Number: 114517

## Product name

Rad51D antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Rad51D Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in COLO 320 cells, Jurkat cells. Observed molecular weight by Western-blot: 23-25kd

## Preparation

This antibody was obtained by immunization of Rad51D recombinant protein (Accession Number: NM\_002878). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:114517(Rad51D antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114517(Rad51D antibody) at dilution of 1:100 (under 10x lens)