

KAZALD1 antibody

Catalog Number: 111917

Product name

KAZALD1 antibody

Specificity

Human; other species not tested.

Antibody description

KAZALD1 Rabbit Polyclonal antibody. Positive IHC detected in human small intestine tissue, human skeletal muscle tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of KAZALD1 recombinant protein (Accession Number: XM_017016715). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111917(KAZALD1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:111917(KAZALD1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:111917(KAZALD1 Antibody) at dilution of 1:50 (under 40x lens)