

ANLN Polyclonal Antibody



Catalog Number: 164917

Product name

ANLN Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

Specificity

Human, Rat

Antibody description

Polyclonal antibody to ANLN

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 845-1124 of human ANLN (NP_061155.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC, IF


Dilutions

Validations



Western blot - ANLN Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ANLN antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Basic Kit .Exposure time: 90s.

Immunofluorescence - ANLN Polyclonal Antibody

Immunofluorescence analysis of A549 cells using ANLN antibody .