

# Anti-Bcl10 Rabbit antibody



Catalog Number: 167917

## Product name

Anti-Bcl10 Rabbit antibody

## Specificity

Human

## Antibody description

Bcl10 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Bcl10

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, ICC/IF, IP


## Dilutions

WB: 1/1000

ICC/IF: 1/50

IP: 1/20

## Validations



Western blot detection of Bcl10 in A549, MCF-7, 293T cell lysates using Bcl10 Rabbit mAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 26kDa.