

LHX4 antibody



Catalog Number: 112217

Product name

LHX4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LHX4 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Observed molecular weight by Western-blot: 43-47 kDa

Preparation

This antibody was obtained by immunization of LHX4 recombinant protein (Accession Number: NM_033343). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112217(LHX4 Antibody) at dilution of 1:300