

# COMMD1 Polyclonal Antibody



Catalog Number: 165217

## Product name

COMMD1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to COMMD1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human COMMD1 (NP\_689729.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions


IF 1:50 - 1:100

## Validations



Immunofluorescence - COMMD1 Polyclonal Antibody

Immunofluorescence analysis of L929 cells using COMMD1 antibody at dilution of 1:100. Blue:



DAPI for nuclear staining.

Immunofluorescence - COMMD1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using COMMD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.