

Sestrin 2 antibody



Catalog Number: 115217

Product name

Sestrin 2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Sestrin 2 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in mouse thymus tissue, HeLa cells, K-562 cells, mouse lung tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 60-66 kDa

Preparation

This antibody was obtained by immunization of Sestrin 2 recombinant protein (Accession Number: NM_031459). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:115217 (SESN2 antibody) at dilution of 1:500



Immunofluorescent analysis of Hela cells, using SESN2 antibody Catalog No:115217 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).