

ASBT antibody

Catalog Number: 108217

Product name

ASBT antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ASBT Rabbit Polyclonal antibody. Positive IP detected in mouse kidney tissue. Observed molecular weight by Western-blot: 38-40 kDa

Preparation

This antibody was obtained by immunization of ASBT recombinant protein (Accession Number: NM_000452). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-SLC10A2 (IP:Catalog No:108217, 4ug; Detection:Catalog No:108217 1:300) with mouse kidney tissue lysate 4000ug.