

Human LILRB1/ILT2/LIR1/MIR7 (His tag) recombinant protein



Catalog Number: 519217

General Information

Gene Name Synonym

Leukocyte Immunoglobulin-Like Receptor Subfamily B Member 1, LIR-1, Leukocyte Immunoglobulin-Like Receptor 1, CD85 Antigen-Like Family Member J, Immunoglobulin-Like Transcript 2, ILT-2, Monocyte/Macrophage Immunoglobulin-Like Receptor 7, MIR-7, CD85j, LILRB1, ILT2, LIR1, MIR7

Protein Construction

A DNA sequence encoding the Human LILRB1/ILT2/LIR1/MIR7 Gly24-His458 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted

samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2 \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\mu\text{g/ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human LILRB1/ILT2/LIR1/MIR7 (His tag)
recombinant protein