

ACE Polyclonal Antibody



Catalog Number: 166617

Product name

Rabbit IgG

ACE Polyclonal Antibody

Applications

Specificity

WB

Human, Mouse

Dilutions

Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to ACE

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-270 of human ACE (NP_690043.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations



Western blot - ACE Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ACE antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.