

General Information

Gene Name Synonym

FIL1 delta; IL-1-related protein 3; Interleukin-1 HY1; Interleukin-1 delta; Interleukin-1 family member 5; Interleukin-1 receptor antagonist homolog 1; Interleukin-1-like protein 1

Protein Construction

A DNA sequence encoding the mature form of human IL1F5 (Q9UBH0) (Met1-Asp155) was expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 99 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human IL1F5 consists of 156 amino acids and predicts a molecular mass of 17.1 KDa. It migrates as an approximately 18 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL1F5 / IL36RN Protein SDS-PAGE