

Anti-HTR1A antibody



Catalog Number: 176317

Product name

Anti-HTR1A antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P, IF

Antibody description

Rabbit polyclonal antibody to HTR1A

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human 5-HT_{1A}. aa range:291-340

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal




Fig1.; Western Blot analysis of various cells using SR-1A Polyclonal Antibody diluted at 1: 500