

TRIM44 antibody



Catalog Number: 116317

Product name

TRIM44 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

TRIM44 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, human liver tissue. Positive IHC detected in human gliomas. Observed molecular weight by Western-blot: 50-55kd

Preparation

This antibody was obtained by immunization of TRIM44 recombinant protein (Accession Number: NM_017583). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116317(TRIM44 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:116317(TRIM44 antibody) at dilution of 1:50 (under 10x lens)