

# H16N3 HA (His Tag) recombinant protein



Catalog Number: 501717

## General Information

### Protein Construction

A DNA sequence encoding the Influenza A virus (A/black-headed gull/Sweden/5/99 (H16N3)) hemagglutinin (AAV91217.1) (Met1-Lys529), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

### Organism

H16N3

### Expression Host

Baculovirus-Insect Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE.

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Asp 20

### Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/black-headed gull/Sweden/5/99 (H16N3)) consists of 521 amino acids and predicts a

molecular mass of 58.6 kDa.

### Formulation

- Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 %glycerol, pH 7.4.
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
  2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Influenza A H16N3 (A/black-headed gull/Sweden/5/99) Hemagglutinin / HA Protein (His Tag)