

General Information

Gene Name Synonym

DCUN1 domain-containing protein 2; Defective in cullin neddylation protein 1-like protein 2

Protein Construction

A DNA sequence encoding the human DCUN1D2 (NP_001014305.1) (Met1-Phe259) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human DCUN1D2 consists of

277 amino acids and predicts a molecular mass of 32.4 KDa. It migrates as an approximately 33 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, 20 % glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide

Storage


Store it under sterile conditions at -20°C to -80°C.

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human DCUN1D2 Protein (His Tag) SDS-PAGE