

RIPK5 antibody



Catalog Number: 114717

Product name

RIPK5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RIPK5 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of RIPK5 recombinant protein (Accession Number: BC060870). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of HepG2 cells, using DSTYK antibody Catalog No:114717 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).