

KDM4A antibody



Catalog Number: 112027

Product name

KDM4A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

KDM4A Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, MCF-7 cells. Observed molecular weight by Western-blot: 55 kDa

Preparation

This antibody was obtained by immunization of KDM4A recombinant protein (Accession Number: XM_005271355). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112027(KDM4A Antibody) at dilution of 1:300



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112027(KDM4A Antibody) at dilution of 1:300