

PMP2 antibody

Catalog Number: 114027

Product name

PMP2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PMP2 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in mouse testis tissue, human brain tissue. Observed molecular weight by Western-blot: 15-20 kDa

Preparation

This antibody was obtained by immunization of PMP2 recombinant protein (Accession Number: NM_002677). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114027(PMP2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114027(PMP2 antibody) at dilution of 1:100 (under 10x lens)