

TEAD1 Polyclonal Antibody



Catalog Number: 165027

Product name

TEAD1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to TEAD1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 135-215 of human TEAD1 (NP_068780.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF, ChIP


Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

ChIP 1:50 - 1:200

Validations



Western blot - TEAD1 Polyclonal Antibody