

# TGIF2LX antibody



Catalog Number: 116027

## Product name

TGIF2LX antibody

## Specificity

human; other species not tested.

## Antibody description

TGIF2LX Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 34kd

## Preparation

This antibody was obtained by immunization of TGIF2LX recombinant protein (Accession Number: BC029920). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:116027(TGIF2LX antibody) at dilution of 1:600