

ZNF74 antibody

Catalog Number: 117027

Product name

ZNF74 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

ZNF74 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive WB detected in Jurkat cells, human brain tissue. Observed molecular weight by Western-blot: 64-72kD

Preparation

This antibody was obtained by immunization of ZNF74 recombinant protein (Accession Number: NM_001256524). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:117027(ZNF74 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells, using ZNF74 antibody Catalog No:117027 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).