

OS9 antibody

Catalog Number: 113427

Product name

OS9 antibody

Specificity

human; other species not tested.

Antibody description

OS9 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 73kd

Preparation

This antibody was obtained by immunization of OS9 recombinant protein (Accession Number: NM_001017956). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113427(OS9 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:113427(OS9 antibody) at dilution of 1:50 (under 10x lens)