

Anti-ZKSCAN4 antibody



Catalog Number: 176427

Product name

Anti-ZKSCAN4 antibody

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ZKSCAN4

Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human zkscan4.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

Validations




Fig1:: Western Blot analysis of MCF7, HepG2, 293, mouse brain, mouse lung cells using ZKSCAN4 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).




Fig2:: Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100