

# TSFM antibody

Catalog Number: 116427

## Product name

TSFM antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TSFM Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human stomach cancer tissue. Positive WB detected in HEK-293 cells, HeLa cells, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of TSFM recombinant protein (Accession Number: NM\_001172696). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116427(TSFM antibody) at dilution of 1:50



(under 10x lens)

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116427(TSFM Antibody) at dilution of 1:600