

Product name

Polyclonal

GAK antibody

Ig Type**Specificity**

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications**Antibody description**

ELISA, WB

GAK Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 150 kDa

Dilutions

Recommended Dilution:


Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of GAK recombinant protein (Accession Number: BC008668). Purification method: Antigen affinity purified.

Validations**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110827(GAK antibody) at dilution of 1:500

Clonality