

RETN Polyclonal Antibody



Catalog Number: 163827

Product name

RETN Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to RETN

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human RETN (NP_001180303.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC


Dilutions

Validations



Immunohistochemistry - RETN Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using RETN antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - RETN Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using RETN antibody at dilution of 1:100 (40x lens).