

# Human HAMP/HEPC/LEAP1/UNQ487/PRO1003 (GST tag) recombinant protein



Catalog Number: 518827

## General Information

### Gene Name Synonym

Liver-expressed antimicrobial peptide 1, Putative liver tumor regressor, HEPC, LEAP1

### Protein Construction

A DNA sequence encoding the Human HAMP/HEPC/LEAP1/UNQ487/PRO1003 Asp60-Thr84 is expressed with a GST tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of PBS□50%glycerol□pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human HAMP/HEPC/LEAP1/UNQ487/PRO1003 (GST tag) recombinant protein