

CAMK2D Polyclonal Antibody



Catalog Number: 160127

Product name

CAMK2D Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CAMK2D

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human CAMK2D (NP_742125.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100

IF 1:20 - 1:50

Validations



Western blot - CAMK2D Polyclonal Antibody

Western blot analysis of extracts of mouse heart, using CAMK2D antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.