

# Anti-CD40 antibody



Catalog Number: 141127

## Product name

Anti-CD40 antibody

## Immunogen

[Rat CD40/TNFRSF5 \(His Tag\) recombinant protein](#)

## Specificity

Rat CD40/TNFRSF5

## Antibody description

Mouse Monoclonal to Rat CD40

## Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rat CD40/TNFRSF5 and conjugated with APC under optimum conditions, the unreacted APC was removed.

## Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

## Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1


## Applications

FCM

## Dilutions

FCM

## Validations



Rat CD40/TNFRSF5 Flow Cytometry (FC) 15893