

Anti-CD40 antibody



Catalog Number: 141127

Product name

Anti-CD40 antibody

Immunogen

[Rat CD40/TNFRSF5 \(His Tag\) recombinant protein](#)

Specificity

Rat CD40/TNFRSF5

Antibody description

Mouse Monoclonal to Rat CD40

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Rat CD40/TNFRSF5 and conjugated with APC under optimum conditions, the unreacted APC was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze !

Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Rat CD40/TNFRSF5 Flow Cytometry (FC) 15893