

CD27 antibody

Catalog Number: 107127

Product name

CD27 antibody

Specificity

Human; other species not tested.

Antibody description

CD27 Mouse Monoclonal antibody. Positive WB detected in Ramos cells, Daudi cells. Positive IHC detected in human tonsillitis tissue. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of CD27 recombinant protein (Accession Number: NM_001242). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


Validations



Ramos cells were subjected to SDS PAGE followed by western blot with Catalog No:107127(CD27 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107127(CD27 Antibody) at dilution of 1:200 (under 10x lens). Heat-mediated antigen retrieval was performed with Tris-EDTA buffer (pH 9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:107127(CD27 Antibody) at dilution of 1:200 (under 40x lens). Heat-mediated antigen retrieval was performed with Tris-EDTA buffer (pH 9).

CD27 antibody



Catalog Number: 107127

with Tris-EDTA buffer(pH9).