

# Enapotamab reference antibody



Catalog Number: 177127

## Product name

Enapotamab reference antibody

## Specificity

Human AXL / UFO

## Antibody description

Anti-AXL / UFO Antibody (Enapotamab)

## Preparation

Recombinant expression and purified from CHO cells.

## Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

## Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal


## Ig Type

IgG1

## Applications


ELISA

## Validations




## SDS-PAGE

Anti-AXL / UFO Antibody (enapotamab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



## SEC-HPLC

The purity of Anti-AXL / UFO Antibody (enapotamab) is 99.31%, determined by SEC-



## HPLC

Bioactivity: ELISA

Immobilized human Axl FC at 2 ug/mL can bind Anti-AXL / UFO Antibody (enapotamab), EC50=0.02735 ug/mL.