

Anti-Prolactin/PRL Mouse antibody



Catalog Number: 169127

Product name

Anti-Prolactin/PRL Mouse antibody

Specificity

Human

Antibody description

Prolactin/PRL Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human PRL (AA: 65-173) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b


Applications

IHC-P, ICC/IF, FC

Validations




Immunofluorescence analysis of HeLa cells using



PRL mouse mAb (green). Blue

Western blot analysis using PRL mAb against




human PRL (AA


Flow cytometric analysis of HeLa cells using PRL mouse mAb (green) and negative control (red).

Anti-Prolactin/PRL Mouse antibody


Catalog Number: 169127



Western blot analysis using PRL mAb against HEK293 (1) and PRL (AA



Immunohistochemical analysis of paraffin-embedded bladder cancer tissues using PRL mouse mAb with DAB staining.



Immunohistochemical analysis of paraffin-embedded muscle tissues using PRL mouse mAb with DAB staining.