

UBE2H antibody



Catalog Number: 116527

Product name

UBE2H antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UBE2H Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse brain tissue, mouse kidney tissue. Positive IHC detected in human colon tissue. Observed molecular weight by Western-blot: 21kd

Validations




Preparation

This antibody was obtained by immunization of UBE2H recombinant protein (Accession Number: NM_003344). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116527(UBE2H antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Immunohistochemistry of paraffin-embedded human colon slide using Catalog No:116527(UBE2H Antibody) at dilution of 1:50

Dilutions

Recommended Dilution:

WB: 1:500-1:5000



Immunohistochemistry of paraffin-embedded human colon slide using Catalog No:116527(UBE2H Antibody) at dilution of 1:50