

# Anti-Annexin VI Rabbit antibody



Catalog Number: 167527

## Product name

Rabbit IgG

Anti-Annexin VI Rabbit antibody

## Applications

## Specificity

WB, IHC-P, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

Annexin VI Rabbit polyclonal antibody

IHC: 1/200

## Preparation


IP: 1/20

Antigen: Recombinant protein of human Annexin-6

## Validations

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

Western blot detection of Annexin-6 in 3T3,C6 cell lysates using Annexin-6 Rabbit pAb(1:1000 diluted).Predicted band size:76kDa.Observed band size:70kDa.

## Ig Type